
FRESH IDEAS IN DIGITAL IMAGING



THE BEST OF BOTH WORLDS







• Small Tabletop Design

- High Resolution Scans in a Fraction of the Time
- Automatic In-line Eraser
- Powerful, Easy-To-Use Software Included
- Value-Priced for any Dental Practice





Apixia Inc. Company Address: 18688 San Jose Ave. Industry, CA 91748 USA Office: +1-626-810-9201 Fax: +1-626-810-9207 Email: info@apixia.com www.apixia.net

Apixia Japan Company Address: 1324, Ooida, Daian-cho, Inabe-shi, Mie-ken 511-0274, Japan Office: +81-59-441-3080 Fax: +81-59-441-3081 Email: sales3@apixia.asia www.apixia.co.jp



TECHNICAL SPECIFICATIONS

COMPUTER REQUIREMENTS*

Windows 11/Windows 10/Windows 7 (32-bit, 64-bit)

USB 2.0/USB 3.0

CD/DVD Drive

*Computer not included

DIGIREX® PSP SCANNER SPECIFICATIONS

Dimensions	183 mm x 138 mm x 258 mm
Weight	3.6 kg
Voltage	AC 100-240V, 50-60Hz
Resolution	17 lp/mm
Connection Type	USB 2.0
Plate Sizes	Size #0, #1, #2, #3
Plate Erasing	Automatic
Software	TWAIN Compatible / DICOM 3.0

INSTANT





LASER FAST IMAGE PROCESSING

Imaging is fast and simple. The smart top-loading design harnesses the power of gravity for simple loading and fuss-free use. In just a short time, the crystal clear radiograph appears on the computer screen, ready for viewing. Simply insert the next exposed plate when ready.

INTELLIGENT



FLEXIBLE IN MORE WAYS THAN ONE

The DIGIREX® PSP Scanner packs tremendous imaging power into a small package. Tiny enough to fit into the most crowded dental operatories for handy chairside use, yet robust enough to be placed in a central location for use by the entire office, the DIGIREX® PSP Scanner is the ideal product for any size dental practice. Imaging with Phosphor plates offers all the flexibility and patient comfort of film, combined with the immediacy and environmental benefits of digital radiography.

COMFORTABLE

REUSABLE

PHOSPHOR PLATE

4 optional plate sizes allows for optimum size selection for each patient



INTUITIVE

POWERFUL SOFTWARE INCLUDED

Apixia PSP Digital Imaging Software is included at no charge with the DIGIREX® PSP Scanner. Packed with value, Apixia includes unlimited user licenses to install on as many computers as you wish, free technical support for as long as you own the product, and complimentary software upgrades for life. Simplicity, ease of use, and an intuitive design were just a few of the things customers told us they wanted in their imaging software.



The software delivers dozens of powerful features,

presented with simplicity and technological elegance that

make the system a snap for every user.

INTERFACE

INTEGRATION

An optional bridge module seamlessly integrates all popular brands of Practice Management Software to Apixia Digital Imaging Software, eliminating the need for double entry of patient data. Just one click and you'll be ready to take radiographs!

Free networking solution included at no additional cost. The Apixia networking software gives you the ability to access all your patient files from any computer in your

*Optional 3rd party practice management bridging may be available in your area. Contact your local Apixia dealer for

